

GrainMapper3D Application Note

Mapping many facets of quartz

Quartz and Dauphiné twinning

Quartz (SiO_2) is one of the most abundant minerals in Earth's crust and is found in nearly all rocks of igneous, metamorphic, and sedimentary origin. It also precipitates from hydrothermal fluids passing through open spaces in subsurface rocks, allowing formation of well-developed crystals like those in **Figure 1**.

Quartz is one of many minerals that develops twinning, or symmetrical intergrowths. One of its most common twin types grows according to the Dauphiné Law but is rarely observed by the unaided eye. The Dauphiné twin has impacted quartz's technological role in society for nearly a century because it disrupts the mineral's piezoelectricity, a property that makes quartz an ideal resonator in oscillator circuits. Quartz is found in everyday electronics, such as clocks, radios, computers, and cellphones. In the oscillator-plate industry, Dauphiné twin removal has been intensely investigated for engineering purposes. In contrast, twins have become a useful tool to reconstruct geological processes.

Dauphiné twins in 3D

Until now, studies of Dauphiné twins have required destructive analysis like acid etching or EBSD to see them, restricting observations of their complex morphologies to two dimensions. **Figure 1** shows a cluster of hydrothermal quartz crystals that each contain a Dauphiné twin, which was measured and characterized in 3D by lab-based DCT [1]. Indexing crystal faces, whose hkl are determined from DCT data, verifies the Dauphiné Law, which transposes positive r $\{10\bar{1}1\}$ and negative z $\{01\bar{1}1\}$ rhombohedral forms. This is seen in **Figure 1** as a "swap" of r and z planes across the curved twin boundary where it intersects an external rhombohedral face. The top and




Figure 1 Correlative maps of quartz obtained by combining lab-based diffraction contrast tomography (DCT) and absorption contrast tomography (ACT). From top to bottom: Crystal facets (gray) are observed by ACT. The hkl families are then assigned to crystal faces, which belong to rhombohedral (r , z) and prism (m) forms. The IPF map shows crystal orientations; Dauphiné twin domains are separated by curved boundaries in each crystal. Distributions of fluid (blue) and solid (black) inclusions are viewed with respect to twin and grain boundaries.

